

Autoadapt Mobility Solutions in partnership with Toyota



INVENTIONS FOR QUALITY OF LIFE

Autoadapt mobility solutions in partnership with Toyota



From June 2006, Autoadapt has entered into a partnership with Toyota for the introduction of their Mobility Programme.

The Toyota Mobility Programme consists of two specific offers:
on one side, the Toyota genuine passenger swivel seat and the passenger lift-up seat
on the other side, adaptations are offered through partnership with Autoadapt and the selected Autoadapt dealers.

Combining Toyota genuine mobility seats alongside adaptation products manufactured by Autoadapt Sweden, the Toyota Mobility Programme offers reliable solutions to meet individual needs.

All products offered by Autoadapt for the Toyota Mobility Programme are tested and approved according to ISO standards and meet current safety regulations for equipment installed in vehicles. Autoadapt is responsible for the installation of the adaptation products and has therefore selected a network of Autoadapt dealers that are geographically spread and have deep knowledge and experience in installing our products.

Autoadapt in partnership with Toyota

Through a partnership with Toyota, Autoadapt can provide Toyota customers with a complementary range of mobility products and adaptations not available at Authorised Toyota Centres.

Autoadapt is a Swedish manufacturer of mobility adaptation products which focuses on quality from design, through engineering, product development and customer service.

The Autoadapt Pan-European network, selected carefully by Autoadapt Sweden, is working with Toyota Centres to propose quality solutions. Toyota mobility centres can guide customers in the logistics of obtaining the adaptations fitted by Autoadapt. Customers can contact their Authorised Toyota Centre for a list of the nearest Autoadapt centres. The Autoadapt products are covered by a 2 years warranty.

Carolift

An electrically powered wheelchair hoist operated with a handheld unit for lifting and lowering the wheelchair as well as for swivelling it in and out of the vehicle.



Roof Spider

A storage box on your rooftop to transport a wheelchair. The Roof Spider can be loaded from either side of the vehicle so it makes no difference if you are a driver or a passenger. Just press the button that corresponds to the side of the Roof Spider that you want to access by a wired remote control.

The design is ideal for foldable, manual wheelchairs up to 25 kg. By storing the wheelchair on the roof you save the space in the back for luggage, groceries etc.



Carony

The Carony gives the disabled passenger the facility to get in and out of the vehicle comfortably.

Dock the wheelchair to the swivel base in the vehicle. Release the seat and slide it over to the swivel base. Now you have made a comfortable transfer without lifting the passenger.



The Carony models



Carony 24" is equipped with easily removable 24" rear wheels. For people who want to move about independently indoors as well as outdoors.

Carony Kids is specially designed to fit a child. The 40 cm deep seat can be adjusted as the child grows. It is available in the 12" or the 24" rear wheel versions.

Carony 12" is a manual wheelchair, equipped with 12" rear wheels. It is a transportation wheelchair, suitable for people who have assistance all the way.

Carony Fixed 12" is a manual wheelchair specially designed to be combined with Turny in high-floored vehicles.



Carony Go

Carony Go is a powered wheelchair. The only crash tested, ISO proofed and CE-labelled powered wheelchair that enables you to travel in a standard small car. It is approved by the Swedish Handicap Institute. The manoeuvre system is produced by Dynamic.

The Carony system is only applicable with the Autoadapt swivel seat Turnout or the seat lift Turny

To be able to use a Carony (page 5) you need to have a Turnout swivel or a Turny mounted in your vehicle. The Carony seat is used both as a wheelchair seat and as a car seat. This means that the wheelchair occupant, always seated, slides over via rails into the vehicle.

The Turnout and Turny have been crash tested, CE-marked and fulfil ISO standards.

Turnout

Turnout is a universal swivel base for the Carony seat. It is available as manually operated or as powered by using a hand held control.



Turny Orbit

The Carony is used in combination with a Turny, if the vehicle is high-floored. Turny swivels the seat outwards and lowers it outside the vehicle. Turny is universal and can be fitted in any high-floored vehicle.



Use the hand held control to operate Turny. The movement stops as soon as the control button is released.



Products with focus on the driver

The products are EMC tested and CE-marked



The "K 5" Gas Ring

This unique Gas Ring offers five different ways of operating the accelerator.

Either:

- (1) push the ring with the back of the fingers or
- (2) push the ring with the thumb from the inside or
- (3) push the ring with the thumb from the outside while holding the hand on the brake lever or
- (4) pull the ring toward the steering wheel with the fingers or
- (5) If you place the special optional "K 5 Paddle Pad" on the steering wheel you can push this to accelerate the car.

Brake Lever

The Brake Lever is mounted by the steering wheel. You brake the car by pushing down the lever and you won't need to take your hand off the steering wheel. It works smoothly with a short stroke and all strength used is multiplied when transmitted to the brake pedal. The handle has integrated buttons for the Hill Holder and the horn.

Remote Transmitter

The steering wheel mounted Remote Transmitter can, when not in use, easily be removed leaving only the holder on the steering wheel. The transmission works via radio waves and activates the electrical functions: position lights, low beam lights, high beam lights, direction indicators, windscreen wipers (two speeds)

Carospeed

Floor mounted hand control unit for accelerator and brake, through a mechanical link. Suitable for people with impaired leg function. A set of primary control functions can be integrated into the unit: indicators, hazard lights, wind-screen washer and wiper, head-lights, etc.

Steering device

Suitable for drivers who steer with one hand and use a hand control unit for acceleration and braking. There are a variety of handgrip alternatives available for this device.

Joey Lift

The Joey platform lift demands very little effort from the user - just drive on to the platform, secure the wheelchair or scooter, press the hand control button and the job is done.

This is the ideal wheelchair stower for people with reduced muscle abilities. The Joey lift will make it very easy for you to bring the wheelchair/scooter along. The Joey has a capacity of loading wheelchairs or scooters weighing up to 159 kg.



Care belt system Careva

Careva stabilisation belts are very helpful for people who need extra support when sitting, for example if the upper body leans to one side or the lower body slides forward in the chair.

The belt system comprises a spinal belt which is fastened vertically over the backrest. The shoulder harnesses, thoracic belt, pelvic belt or a combi harness are then installed in the grooves of the spinal belt.



Toyota Genuine Mobility Programme

The Toyota Genuine Mobility Programme consists currently of a passenger swivel seat and a passenger lift-up seat:

Swivel seat

The Toyota swivel seat provides a comfortable and easy way to get in and out of vehicles. A swivel base is combined with either a low construction seat or the modified original seat. The seat can be swivelled outwards manually over the doorsill with a lever that is attached to the base. For lower roof height vehicles, the swivel base is attached to a specially designed low construction seat, which is a thinner version of the standard seat.

If headroom is not a concern, the swivel base is attached to the original seat, which only needs to be slightly modified. In either case, mobility adaptations will not alter the comfort or the interior style of your vehicle.



Lift-up seat

The Toyota genuine lift-up seat makes getting in and out of vehicles considerably easier. The lift-up seat is a powered, outwards swivelling car seat lift, specially designed for Toyota vehicles that have a higher step up. The lift-up seat comes with the same fabric as the original standard seat to match the vehicle's interior. The lowering and lifting of the seat is operated with a switch on the side of the seat or with a wired remote control. The sliding and the reclining of the seat inside the vehicle are also operated with switches situated on the side of the seat.



The Toyota Network

Both mobility seats are ordered and installed at your Authorised Toyota Centre. Original seat anchorages are used to install the mobility seats in Toyota Vehicles. All genuine Toyota seats are thoroughly tested according to applicable ECE directives, match with the OE fabric and come with a three-year warranty.

For further details and to test the Toyota genuine mobility products, contact your local Authorised Toyota Centre.

The Autoadapt selected dealer network



GLENEAGLES CONVERSIONS

95 Bonnygate
Cupar
Fife, KY15 4LG
Tel. 01334 657 722
Fax. 01334 657 711
sales@gleneaglesconversions.co.uk
8am - 5pm

Mc ELMEEL MOBILITY SERVICES

15 Ballyscandal Road
Armagh, BT61 8BL
Tel. (028) 3752 5333
Fax. (028) 3752 7740
info@mobility-services.com
9am - 5.30pm Mon-Thur
9am - 4.30pm Fri

MOBILITY CONVERSIONS

Langley Road South
Salford
Manchester, M6 6SY
Tel. 0161 745 9737
Fax. 0161 743 0190
info@mobilityconversions.co.uk
9am - 5pm

CITY MOTOR SERVICES

Unit C10
Springmeadow Business Park
Springmeadow Road
Cardiff
South Glamorgan, CF3 2ES
Tel. 02920 364 800
Fax. 02920 364 900
info@citymotorservices.co.uk
9am - 5.30pm

ADAPTACAR

Pathfields Business Park
South Molton
Devon, EX36 3LH
Tel. 01769 572 785
Fax. 01769 574 984
info@adaptacar.co.uk
8.30 - 5pm

SVO

Redman Road
Porte Marsh Ind Estate
Calne
Wiltshire, SN11 9PR
Tel. 01249 815141
Fax. 01249 822241
info@specialisedvehicleoptions.com
8.30 - 5pm

PB CONVERSIONS

Unit 6, Acacia Close
Cherrycourt Way
Leighton Buzzard
Bedfordshire, LU7 4QE
Tel: 01525 850588
Fax: 01525 850588
info@pbconversions.co.uk

STEERING DEVELOPMENTS

Unit, 5 Eastman Way
Hemel Hempstead
Herts, HP2 7HF
Tel. 01442 212 918
Fax. 014442 240254
enquiries@steeringdevelopments.co.uk
8.30am - 5pm Mon-Thur
8.30am - 4pm Fri

KEWAL MOBILITY

Scotchburn Garth
Skerne Road
Driffield
East Yorkshire, YO25 6EF
Tel. 0845 6446004
Fax. 01377 252500
keith@kewal.com
9am - 5pm

DES GOSLING MOBILITY

Unit 10, Station Yard
Station Road
Melbourne
Derbyshire, DE73 1BQ
Tel. 01332 863 742
Fax. 01332 863 742
info@desgosling.wanadoo.co.uk
8.30am - 5pm

TECH MOBILITY ESSEX

Purley Farm Barns
Off Colne Road
Coggeshall
Essex
CO6 1TH
Tel. 01376 564499
Fax. 01376 564992
phil@techmobilityuk.com
8am - 5pm

TECH MOBILITY SUSSEX

Unit 1, Crosspost Industrial Park
Cowfold, Bolney
West Sussex, RH17 5QU
Tel. 01444 882233
Fax. 01444 882800
sales@techmobility.co.uk
8.30am. - 5.30pm





Autoadapt AB is a Swedish Car Adaptation Company with more than 20 years of experience in car adaptations, - designing, producing and distributing products for Car Adaptation world wide. Products, such as the Carony car seat docking system, in combination with the Turnout swivel base, or the Turny car seat lift have made the company well known in our industry.

Autoadapt's main ambition is to develop well designed products that are of the highest quality, safe, CE-marked and tested according to current laws and regulations.



Crash test of a Turny Orbit in a Toyota.

For more
information
www.autoadapt.co.uk



Autoadapt UK Ltd

Unit 4, Windsor Industrial Estate
Rupert Street, Aston
Birmingham B7 4PR

Phone: 0121 333 5170
E mail: contact@autoadapt.co.uk
www.autoadapt.co.uk
www.autoadapt.com